

Product datasheet for **RK300057**

MIR146b Rat qPCR Template Standard (MI0006342)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGAGAACTGAATTCCATAGGC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0005595

Synonyms: MI0006342; MIMAT0005595; MIR146b rno-mir-146b; rno-miR-146b-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)