

## Product datasheet for **RK300072**

### MIR181b-2 Rat qPCR Template Standard (MI0000927)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AACATTCATTGCTGTCGGTGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000859
Synonyms:	MI0000927; MIMAT0000859; MIR181b-2 rno-mir-181b-2; rno-miR-181b-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)