

## Product datasheet for **RK300096**

### MIR196c Rat qPCR Template Standard (MI0006136)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGGTAGTTTCGTGTTGTTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0005303
Synonyms:	MI0006136; MIMAT0005303; MIR196c rno-mir-196c; rno-miR-196c-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)