

Product datasheet for **RK300100**

MIR19b-2 Rat qPCR Template Standard (MI0000848)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGTTTTGCAGATTGCGTTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017097
Synonyms:	MI0000848; MIMAT0017097; MIR19b-2 rno-mir-19b-2; rno-miR-19b-2-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)