

Product datasheet for **RK300126**

MIR218a-2 Rat qPCR Template Standard (MI0000957)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TTGTGCTTGATCTAACCATG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000888

Synonyms: MI0000957; MIMAT0000888; MIR218a-2 rno-mir-218a-2; rno-miR-218a-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)