

Product datasheet for **RK300128**

MIR219a-2 Rat qPCR Template Standard (MI0000960)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGATTGTCCAAACGCAATTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000889
Synonyms:	MI0000960; MIMAT0000889; MIR219a-2 rno-mir-219a-2; rno-miR-219a-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)