

Product datasheet for **RK300136**

MIR23b Rat qPCR Template Standard (MI0000853)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GGGTTCCCTGGCATGCTGA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017099
Synonyms:	MI0000853; MIMAT0017099; MIR23b rno-mir-23b; rno-miR-23b-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)