

Product datasheet for **RK300141**

MIR26b Rat qPCR Template Standard (MI0000858)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TTCAAGTAATTCAGGATAGGT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000797

Synonyms: MI0000858; MIMAT0000797; MIR26b rno-mir-26b; rno-miR-26b-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)