

Product datasheet for **RK300144**

MIR28 Rat qPCR Template Standard (MI0000861)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | AAGGAGCTCACAGTCTATTG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0000800 |
| Synonyms: | MI0000861; MIMAT0000800; MIR28 rno-mir-28; rno-miR-28-5p |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)