

## Product datasheet for **RK300149**

### **MIR293 Rat qPCR Template Standard (MI0012609)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** ACTCAAACGTGTGACACT

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0012847

**Synonyms:** MI0012609; MIMAT0012847; MIR293 rno-mir-293; rno-miR-293-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)