

## Product datasheet for **RK300156**

### MIR2985 Rat qPCR Template Standard (MI0015469)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGTTATAGTATCCCACCTAC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017891
Synonyms:	MI0015469; MIMAT0017891; MIR2985 rno-mir-2985; rno-miR-2985
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)