

Product datasheet for **RK300157**

MIR299a Rat qPCR Template Standard (MI0000970)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGGTTTACCGTCCCACATAC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000901

Synonyms: MI0000970; MIMAT0000901; MIR299a rno-mir-299a; rno-miR-299a-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)