

Product datasheet for **RK300211**

MIR349 Rat qPCR Template Standard (MI0000636)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CAGCCCTGCTGTCTTAACCT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0000599

Synonyms: MI0000636; MIMAT0000599; MIR349 rno-mir-349; rno-miR-349

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)