

Product datasheet for **RK300224**

MIR3549 Rat qPCR Template Standard (MI0015405)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	ATCTTGCCCGCAAAGACC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017807
Synonyms:	MI0015405; MIMAT0017807; MIR3549 rno-mir-3549; rno-miR-3549
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)