

Product datasheet for **RK300230**

MIR3556b Rat qPCR Template Standard (MI0015414)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACCGATTTCAAATGGTGCTA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017818

Synonyms: MI0015414; MIMAT0017818; MIR3556b rno-mir-3556b; rno-miR-3556b

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)