

## Product datasheet for **RK300234**

### MIR3560 Rat qPCR Template Standard (MI0015422)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** CAAATCCTTGCCCGGTGC

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0017829

**Synonyms:** MI0015422; MIMAT0017829; MIR3560 rno-mir-3560; rno-miR-3560

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)