

Product datasheet for **RK300247**

MIR3576 Rat qPCR Template Standard (MI0015448)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACACACCAAGGATAATTTCTC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017862

Synonyms: MI0015448; MIMAT0017862; MIR3576 rno-mir-3576; rno-miR-3576

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)