

Product datasheet for **RK300251**

MIR3580 Rat qPCR Template Standard (MI0015453)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGCTTGACTTATTCTGGTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017868

Synonyms: MI0015453; MIMAT0017868; MIR3580 rno-mir-3580; rno-miR-3580-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)