

Product datasheet for **RK300257**

MIR3588 Rat qPCR Template Standard (MI0015466)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TCACAAGTTAGGGTCTCAGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017887

Synonyms: MI0015466; MIMAT0017887; MIR3588 rno-mir-3588; rno-miR-3588

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)