

Product datasheet for **RK300267**

MIR3596c Rat qPCR Template Standard (MI0015459)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACTATACAACCTCCTACCT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017877

Synonyms: MI0015459; MIMAT0017877; MIR3596c rno-mir-3596c; rno-miR-3596c

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)