

Product datasheet for **RK300276**

MIR375 Rat qPCR Template Standard (MI0006140)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: GCGACGAGCCCCTCGACA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017301

Synonyms: MI0006140; MIMAT0017301; MIR375 rno-mir-375; rno-miR-375-5p

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)