

Product datasheet for **RK300290**

MIR409b Rat qPCR Template Standard (MI0015454)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGGGGTTACCGAGCAACA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0017870

Synonyms: MI0015454; MIMAT0017870; MIR409b rno-mir-409b; rno-miR-409b

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)