

## Product datasheet for **RK300299**

### **MIR433 Rat qPCR Template Standard (MI0001724)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Forward Sequence:</b>	TACGGTGAGCCTGTCATTA
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0017192
<b>Synonyms:</b>	MI0001724; MIMAT0017192; MIR433 rno-mir-433; rno-miR-433-5p
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
<b>Quality Control:</b>	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)