

## Product datasheet for **RK300314**

### MIR483 Rat qPCR Template Standard (MI0003485)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** AAGACGGGAGAAGAGAAGGG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0017201

**Synonyms:** MI0003485; MIMAT0017201; MIR483 rno-mir-483; rno-miR-483-5p

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)