

Product datasheet for **RK300334**

MIR539 Rat qPCR Template Standard (MI0003526)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | GGAGAAATTATCCTTGGTGTG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0003176 |
| Synonyms: | MI0003526; MIMAT0003176; MIR539 rno-mir-539; rno-miR-539-5p |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)