

## Product datasheet for **RK300354**

### MIR6318 Rat qPCR Template Standard (MI0021840)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CTGCCTGGCGCAGGGCCT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0025055
Synonyms:	MI0021840; MIMAT0025055; MIR6318 rno-mir-6318; rno-miR-6318
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)