

Product datasheet for **RK300361**

MIR6324 Rat qPCR Template Standard (MI0021848)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TCAGTAGGCCAGACAGCAAG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025063

Synonyms: MI0021848; MIMAT0025063; MIR6324 rno-mir-6324; rno-miR-6324

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)