

Product datasheet for **RK300397**

MIR761 Rat qPCR Template Standard (MI0012615)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GCAGCAGGGTGAAACTGACA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0012853
Synonyms:	MI0012615; MIMAT0012853; MIR761 rno-mir-761; rno-miR-761
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)