

## Product datasheet for **RK300401**

### MIR7a-2 Rat qPCR Template Standard (MI0000836)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGAAGACTAGTGATTTTGTT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000606
Synonyms:	MI0000836; MIMAT0000606; MIR7a-2 rno-mir-7a-2; rno-miR-7a-5p
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)