

Product datasheet for **RK300407**

MIR875 Rat qPCR Template Standard (MI0012604)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TATACCTCAGTTTATCAGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0012842

Synonyms: MI0012604; MIMAT0012842; MIR875 rno-mir-875; rno-miR-875

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)