

## Product datasheet for **RN200088**

### RT1-CE5 (NM\_001033986) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** RT1-CE5 (NM\_001033986) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** RT1-CE5  
**Synonyms:** RT1-C/E5; RT1-E2; rt1-EI  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200088 representing NM\_001033986  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGGGCGATGGCACCGCGCACGCTGCTCCTGCTGCTGGCGGCGCCCTGGCCCAACCCGGACCCGCG  
CGGGCTCACACTCGTGCCTATTTCCACACCGCCGTGTCCGGCCCGCCCTCGGGAGCCCGGTTTCAT  
TATCGTCGGCTATGTGGACCACACGGAGTTCGTGCGCTACGACAGCGACGACAGAGAAATCCGAGATACGAG  
CCGCGGGCGCGGTGGATGGAGCGGGAGGGCCGGAGTATTGGGAGCAGCAGACACAGATCGTCAAGGGAC  
ACGAATGGAGTTTCCGAGTGAGCCTGAGGAATCTGCTAGGCTACTACAACCAGAGCGAGGGCGGCTCTCA  
CACACTCCAGGTGATGTCTGGATGTGACATGGGGTCCGACCAGAGCCTGCGCCCGGGTATGAGCAATAC  
GCCTATGACGGCCGCGATTACATCTCCCTGAACGAAGACCTGAAGACGTGGTCCGGCTGCGGACTTGGCAG  
CACAGATCACTCTAAACAAGTGGGATCAGGCTGGTGAATTAGACTATTACAAGGCCTACCTGGAGGGTCC  
ATGCCTGGAGTGGCTTCTCAGATACCTGGAGCACGGGAAGGAGACACTGATGCGCTCAGATCCCCAAAG  
GCACATGTGACCCTTACCCCAGACCTGAAGGTGATGTGACCCTGAGGTGCTGGGCCCTGGGCTTCTACC  
CTGCTGACATCACCTGACCTGGCAGTTGAATGGGGAGGACCTGACCCAGGACATGGAGCTTGTGGAGAG  
CAGGCCTGCAGGGGATGGAACCTTCCAGAAGTGGGCATCTGTGGTGGTGCCTCTTGGGAAGGAGCAGAAT  
TACACATGCCTTGTGGAGCATGAGGGGCTGCCTGAGCCCTCACCTTTAGATGGGAGGCTCCTCCCTCCA  
CCAACCTCAACACGGGAATGTCTGTTATTCTTGGAACTGTGACCATCATTGCAGTTATGGCCATCATTGT  
AGTTGTAGCCGTCATTAGAACTGTTGTGAGGAAGAGATGGATAAAAAACAGCTTTTCTTCTACAAGTGG  
AAAGGAGTCTACACCCCTGCTCCAGGCAGGACAGCTCCAGAGCTCTGATGTGTCTCTCCAGATTGTA  
AAGCATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_001033986
<b>Insert Size:</b>	1128 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001033986.1, NP_001029158.1</u>
<b>RefSeq Size:</b>	1634 bp
<b>RefSeq ORF:</b>	1128 bp
<b>Locus ID:</b>	309607
<b>Cytogenetics:</b>	20p12