

## Product datasheet for **RN200128**

### Yif1b (NM\_001014810) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Yif1b (NM\_001014810) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Yif1b  
**Synonyms:** Yip1b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200128 representing NM\_001014810  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCCACGTCCGGTAGGGGCGCATCGCGCCTCAAAGCGGAGGGTCCCAGTGTCCAACCGGGCATGG  
CTGATCCCCACCAGTTTTTTGACGACACGAGTTCAGCCCCGAGCCGGGGTATGGAGGCCAGCCATCACC  
TGGTAGCCTGGGTACCCACCTCTTCATCGGAGGCAGCCTTCTGGCTGCTCCTATGTCCAACATGGCG  
ATGGCCTACGGGAGCAGCCTGGCTGCTCAGGGCAAAGAGCTGGTAGATAAGAATATCGACCCTTCATCC  
CTGTACAGCAAGCTCAAGTATTACTTCGCCGTGGACACCCTGTACGTGGGCAAGAAGCTCGGACTGCTCGT  
CTTTCCCTACCTTACCAGGACTGGGAGGTCCAGTACCAACAGGACACCCCGTGGCGCCCGCTTTGAC  
ATCAATGTCTCCAGACCTTACATTCAGCCATGGCTTTCATAACCTACATCTTGGTGGCCGGCCTTGAC  
TGGGACCCAGGACAGTTCTCCCAGACCTCTGGGACTACAGGCGAGCTCGGCACTGGCTTGGCTGAC  
TCTGGAAGTCGTAGCCATTCTGCTCAGTCTTTACCTGGTCACTGTCAACACCGACCTTACCACATTGAC  
CTGGTGGCCTTCTGGGTTACAAATATGTGGGATGATTGGCGGGTCTCACAGGCCTGCTCTTTGGGA  
AGATTGGCTACTACCTGGTCTTGCCTGGTGTGTGTCCATCTTTGTGTTTCATGATCCGGACACTGCG  
GCTCAAGATCCTGGCCCAGGACGGCCGAAGGAGTGCCGGTGGCGGGGCTCGGAACCACTGCGCATG  
TACTTGACTATGGCGGTGGCGGCTGCGCAGCCGTAATCATGTACTGGCTCACCTTCCACCTTGTACGGT  
GA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001014810  
**Insert Size:** 912 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001014810.1</a></u> , <u><a href="#">NP_001014810.1</a></u>
<b>RefSeq Size:</b>	1183 bp
<b>RefSeq ORF:</b>	912 bp
<b>Locus ID:</b>	292768
<b>UniProt ID:</b>	<u><a href="#">Q6PEC3</a></u>
<b>Cytogenetics:</b>	1q21
<b>Gene Summary:</b>	Involved in the anterograde traffic pathway from the endoplasmic reticulum to the plasma membrane and the organization of the Golgi architecture (Probable). Plays a key role in targeting to neuronal dendrites receptors such as HTR1A (PubMed:18685031). [UniProtKB/Swiss-Prot Function]