

Product datasheet for **RN200158**

Tp63 (NM_001127341) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tp63 (NM_001127341) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tp63
Synonyms:	Ket; P73l; Tp73l; Trp63
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >RN200158 representing NM_001127341
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATTTTGAACCTTCACGGTGTGCTACCTACAGTACTGCCCTGACCCTTACATCCAGCGTTTCATAG
 AAACCCCATCTCATTCTCCTGGAAGAAAGTTATTACCGGTCCGCCATGTCGCAGAGCACCCAGACAAG
 TGAGTTCTCAGCCAGAGGTGTTCCAGCATATCTGGGATTTTCTGGAACAGCCTATATGCTCAGTACAG
 CCCATCGACTTGAACCTTGTGGACGAACCATCAGAAAATGGTGAACAAACAAGATTGAGATTAGCATGG
 ATTGTATCCGCATGCAAGACTCAGACCTCAGTGACCCCATGTGGCCACAGTACACGAACCTGGGGCTCCT
 GAACGGCATGGACCAGCAGATTCAGAACGGCTCCTCATCTACCAGCCCTATAACACAGACCATGCACAG
 AACAGCGTGACGGCACCCCTCGCCCTATGCACAGCCAGCTCAACCTTCGATGCCCTTTCTCCATCCCCTG
 CCATTCCCTCAAACACAGATTACCCAGGCCACACAGCTTCGATGTGCCTTCCAGCAGTCAAGCACCGC
 CAAGTCAGCTACCTGGAGTATTCCACCGAACTGAAGAACTCTACTGCCAGATTGCAAAGACCTGCCCC
 ATCCAGATCAAGGTGATGACCCACCCACAGGGCGCCGTATTTCGTGCCATGCCTGTCTACAAGAAAG
 CCGAGCATGTACCGAGGTTGTGAAACGATGTCCTAACCCAGAGCTGAGCCGCGAGTTCAATGAGGGACA
 GATTGCCCTCCCAGTCACTGATTTCGAGTAGAAGGGAACAGCCATGCCAGTATGTAGAAGATCCTATC
 ACAGGAAGGCAGAGCGTGCTGGTCCCTTATGAGCCACCACAGGTTGGCACTGAATTCACAACAGTCTGT
 ACAATTTTCATGTGCAACAGCAGCTGTGTCGGAGGAATGAACCGCCGTCCAATTTTAAATCATCGTTACTCT
 GGAAACCAGAGATGGGCAAGTCTGGGCCGACGTTGCTTTGAGGCCCGGATCTGCGCTTGCCAGGAAGA
 GACCGGAAGGCCGATGAAGACAGCATCAGAAAGCAGCAAGTATCAGACAGCGCAAAGAACGGCGATGGTA
 CGAAGCGCCCTTCCGTGAGAATACCCACGGAATCCAGATGACTTCCATCAAGAAACGGAGATCCCAGA
 TGATGAGCTGCTGTACCTACCAGTGAGAGGCCGTGAGACTTATGAAATGCTGCTCAAGATCAAGGAGTCG
 CTCGAGCTCATGCAGTATCTCCCTCAGCACAGTACGAGAGCTACAGGCAGCAGCAGCAGCAGCAGCACC
 AACACCTACTTCAGAAACATCTCCTTTCAGCCTGCTTCAGGAATGAGCTTGTGGAGTCCCGGAGAGAAGC
 TCCGACACAGTCTGACGTCTTCTTAGACATTCCAACCCCAACCCTCAGTGTACCCA**TAG**

ACCGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001127341

Insert Size: 1464 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127341.1](#), [NP_001120813.1](#)

RefSeq Size: 1611 bp

RefSeq ORF: 1464 bp

Locus ID: 246334

UniProt ID: [Q9JJP6](#)

Cytogenetics: 11q22

Gene Summary: This gene encodes tumor protein p63, a member of the p53 family of transcription factors involved in cellular responses to stress and development. The family members include p53, p63, and p73, which have high sequence similarity to one another. This similarity allows p63 and p73 to transactivate p53-responsive genes causing cell cycle arrest and apoptosis. The family members can interact with each other in many ways, including direct and indirect protein interactions. This results in mutual regulation of target gene promoters. Both alternative splicing and the use of alternative promoters result in multiple transcript variants encoding different protein isoforms.[provided by RefSeq, Dec 2009]
Transcript Variant: This variant (3) lacks several 3' exons but has an alternate 3' exon, as compared to variant 1. The resulting isoform (c, also known as TA2KETgamma) has a shorter and different C-terminus, as compared to isoform a.