

Product datasheet for **RN200185**

Stau2 (NM_134466) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Stau2 (NM_134466) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Stau2
Synonyms:	MGC93196
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >RN200185 representing NM_134466
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAAATCCCAAAGAGAAAATCCAATGTGTCTGGTAAATGAGTTAGCCCGTTCCATAGCATCCAAC
 CCCAGTATAAACTTCTGAATGAAAGCGGCCTGCTCATTGGAAGATGTTTTCGGTGCAGCTGAGTCTGGG
 TGAGCAGACATGGGAATCCGAAGGGAGCAGTATAAAGAAGGCCAACAGGCTGTTGCAAATAAAGCTTTG
 ACTGAATCTACGCTTCCCAAACCCGTTCAAGAAACCCAAAAGTAATGTCAATAATAACCCAGGCAGTA
 TAACTCCGACTGTGGAAGTGAATGGGCTTGCTATGAAAAGGGGAGAGCCTGCCATCTACAGGCCGTTAGA
 CCCAAAGCCATCCCAAATTATAGAGCTAATTACAACCTCCGGGGCATGTACAACCAGAGATATCATTGC
 CCAATGCCAAGATCTTTATGTTCAAGTGTGACTGTAGGAAATAATGAATCTTTGGCGAAGGGAAGACTC
 GGCAAGCTGCCAGACACAATGCTGCGATGAAAGCTCTTCAAGCCCTACAGAATGAGCCAATCCAGAAAA
 GTCTCCTCAGAATGGTGAATCAGGAAAAGAAATGGATGACGATAAAGATGCAAATAAATCTGAAATAAGC
 TTAGTGTTTTGAAATTGCGCTGAAGAGAAAATATGCCCGTCAGTTTTGAGGTTATTAAGAAAAGTGGACCAC
 CACATATGAAGAGTTTTGTTACCCGGGTGTCAAGTGGGAGAATTTCTGCAGAAGGAGAGGAAAATAGCAA
 AAAGCTCTCCAAGAAGCGTGCAGCAACCACGGTCTTACAGGAGCTTAAGAAAATCCACCTCTTCTGTG
 ATAGAAAAGCCAAAATTTTTTAAAAACGCCCTAAAACAATAATAAAGGCTGGACCCGAGTATGGTC
 AAGGAATGAACCCATTAGCCGCTGGCTCAAATCCAGCAGGCCAGAAAGGAAAAGGAGCCTGATTACGT
 TTTACTTTCAGAACGTGGAATGCCTCGCCGTCGAGAGTTGTAATGCAGGTCAAGGTAGGCAATGAAGTT
 GCGACTGGAACAGGACCTAATAAAAAGATAGCCAAGAAAATGCTGCAGAAGCTATGCTGTTACAGCTCG
 GCTATAAAGCGTCCACTAGTCTTCAGGATCAGCTCGACAAGACCCGGGAAAACAAAGGATGGAGTGGTCC
 AAAGCCTGGATTTCCGGAACCGCAATAACACTCCAAAAGGAATCTTTCATCTTTCTCCTGATGTTTAT
 CAAGAGATGGAAGCCAGCCGCCACAGAGTGACCTCTGGCACAACCTTAGGCTACTTATCACCCAAAGATA
 TGAACCAACCTTCAAGCTCCTTCTCAGTGTAGAGTCTCCCTCACCGACTAGTTCAGCCCGAGTGCCAG
 GGAACCTCTAATGAATGGAACATCTCCTGCTGCTGAAGCCATAGGTTTAAAAGGAAGTTCTCCTACTCCC
 CCTTGTCTTTCAGTACAGCCTTCAAACAACCTGGAATATTTAGCAAGGATTAAGGCTTTCAGGTA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_134466

Insert Size: 1539 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_134466.2](#), [NP_604461.2](#)

RefSeq Size: 2912 bp

RefSeq ORF: 1539 bp

Locus ID: 171500

UniProt ID: [Q68SB1](#)

Cytogenetics: 5q11

Gene Summary: binds RNA; plays a role in targeting mRNA to neuronal dendrites [RGD, Feb 2006]