

## Product datasheet for **RN200196**

### **Cabp1 (NM\_133529) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cabp1 (NM\_133529) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cabp1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200196 representing NM\_133529  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCAACTGCGTCAAGTCGCCACTAAGAAACCTCTCGAAGAGATGCGCCAGGAGGAGAAGACCAGCT  
ACATGGCGGTGCAGACAAGTGAGGATGGGCTGGCGGATGGCGGTGAGCTCCCTGGACCCTCATGATGCT  
GGCCCAGAAGTGTGCGGTATGCACAACCTGCTGGTCCCGCTGCATCTTCTCCGAAAGGTTTTGCC  
GAGAACAGACAACCTGATAGATCGTTACGACCAGAGGAGATCGAAGAGCTCCGAGAGGCCTTCCGGGAAT  
TTGACAAGGACAAAGATGGCTACATCAACTGCCGGACCTGGGGAAGTGCATGCGCACCATTGGGCTACAT  
GCCCACTGAGATGGAAGTGCATCGAGCTGTCCAGCAGATCAACATGAACCTGGGTGGCCAGCTGGATTTC  
GATGACTTTGTGGAAGTGCATGGGCCCTAAACTCCTGGCGGAGACAGCAGATATGATCGGTGTAAGGAAAC  
TGAGAGACGCCTTTCGGGAGTTTGACACCAACGGTGACGGGGAGATAAGCACGAGCGAGCTGCGAGAGGC  
CATGAGGAAGCTCCTGGGTCATCAGGTGGGACACCGGGACATAGAGGAAATTATCCGAGACGTGGACCTC  
AATGGAGATGGACGAGTGGACTTCGAAGAGTTTGTCCGGATGATGTCTCGCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_133529  
**Insert Size:** 684 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_133529.2</a></u> , <u><a href="#">NP_598213.1</a></u>
<b>RefSeq Size:</b>	1125 bp
<b>RefSeq ORF:</b>	684 bp
<b>Locus ID:</b>	171051
<b>UniProt ID:</b>	<u><a href="#">O88751</a></u>
<b>Cytogenetics:</b>	12q16
<b>Gene Summary:</b>	<p>component of calcium-mediated signal transduction; may be involved with the dendritic cytoskeleton [RGD, Feb 2006]</p> <p>Transcript Variant: This variant (1) has an alternate 5' sequence, as compared to variant 3. The encoded isoform 1 has a distinct and shorter N-terminus, as compared to isoform 3.</p>