

Product datasheet for **RN200219**

Cflar (NM_001033864) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cflar (NM_001033864) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Cflar
Synonyms:	Flip
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >RN200219 representing NM_001033864
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCCTGAGCACTGTGTCTGCTGAGGTCATTACCAGGTGGAGGAGTCTCTTGATGAGGATGAGAAGG
 AGATGATGCTCTTCTTGTGTAGAGATGTGACTGAGAACCCTGGCTCCACCCAACGTGAGGACCTCCTGGA
 TTGTTTAAAGTGAGAGAGGCCAGCTCTCCTTCGCTGCCTTAGCTGAGCTGTGTACAGGGTGCGGCGTTT
 GACCTGCTCAAGAGGATTTGAAGACAGATAAAGCAGCGGTGGAGGACCATCTGTGCAGAAGCCCTCACC
 TTGTTTCCGATTATAGGGTGTGCTGATGGAGATTGGGGAGAATTTGAATCAATCTGATGTATCCTCCTT
 AATTTTCCTTACAAAGGATTACACAGGCCGAGGCAAGGTAGCCAAGGACAAGAGTTTCTGGACCTGGT
 ATTGAAGTGGAGAACTGAATTTGATTGGTTCCGATCAGTTGAATTTATTAGAAAAATGCCTGAAGAGCA
 TCCACAGAATAGACTTGAAGACAAAAATCCAGAAGTACACGCAGTCCAGCCAAGAAGCAAGATCAAATAT
 GAATGCTCTCCAGCTTCTCTCCAAAATTGAGTATCAAGGAACATTTATATAACTCAAGGCTCCAGAAT
 GGGCGAAGCAAAGAGCCAAGATTTGTGGAACACCATGGCATTAGAGAAAACCGGTGAAGACATCCATTC
 AGGAATCAGGAATTTTTCTGCCTCCGCACATACATGAAGAGAGCTACAGGATGCAGAGCAAGCCTTAGG
 GATCTGCTTATCATCGACTGTATTGGCAATGACACAGAAATCTTCGAGAGACCTTCACTTCTTGGGC
 TATCGCGTCCAGCCTTCTTGTTCCTCAGTTTCCAGTTCACATGAAATAACCCAGATTGTTCCGCCATTCTCCAATA
 TGACCCAACATCAAGACTATGACAGCTTTGATGTGTCTGGTGAGCCGAGGAGGCTCCCAAAGCATGAT
 GGGCGTAGATCAGGTTTACTCAGGATTCCTTGGAGAATGTCAAGGACATGTTCAAGGGGGACATGTGC
 CCTTCTCTCAGAGGAAAACCCAAGCTCTTTTTTATTCAAACATGAGGGCTTGGAGATAACAGCCTGG
 AAGTAGATGGCCAGCAATAAAAAATGTAACCTAGGCACCTGACATCCAGACTGCACAACATCACC
 GGACGCTGATATCTTCTGGAGCCTGTGCACAGCAGACGTATCCCGTTGGAGCAGCCCTCCAGCTCGCTC
 TCTGTGTACCTGCAGAACTCTCCAGCTGCTGGAGCAGGGCAGGAAGCGCTCACTCGTGGATCTCCATG
 TTGAAGTATGGACAAAGTGTATGCCTGGAACAGTAATGTTCTTCCAAGGAGAAATACTACCTCAGTCT
 GCAGCAGACTCTGAGGAAGAACTCATCTGTCTCACACATGA

ACGGTACGGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001033864
Insert Size: 1443 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033864.2</u> , <u>NP_001029036.1</u>
RefSeq Size:	1799 bp
RefSeq ORF:	1443 bp
Locus ID:	117279
Cytogenetics:	9q31
Gene Summary:	an anti-apoptosis protein [RGD, Feb 2006]