

## Product datasheet for **RN200228**

### Calcr (NM\_053816) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Calcr (NM_053816) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Calcr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >RN200228 representing NM\_053816  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGACTCCACGGAGATCCAGAATGAAAAGCGGAACTCCGAAAACCGAAGATGAGGTTCTTCTCCTGA  
 ACAGGTTACCCCTGCTGCTCCTGCTCCTAGTGAGCCCTACTCCAGTTCTTCAGGCTCCTACCAATCTCAC  
 TGA**CTCCGGCCTTGATCAGGAGCCATTCCTGTACTTGGTTGGACGCAAGAAGCTTCTGGATGCTCAGTAC**  
 AAGTGCTATGACCGGATTCAGCAGTTGCCCCCTTATGAGGGAGAAGGTCCATATTGCAACCGAACTTGGG  
 ATGGATGGATGTGCTGGGACGACTCCAGCTGGTGAATGTCCTATCAGCACTGCCCTGACTACTTTCC  
 GGACTTTGACCCAACAGAAAAGTTTCAAATACTGCGATGAAAATGGAGAGTGGTTTAGACACCCTGAC  
 AGCAACCGAACCTGGTCCA**ACTACTCTGTGCAATGCCTTCACTCCTGACAACTGCATAATGCATACG**  
 TTCTTTATTACCTAGCACTTGTGGGCCACTCCATGTCAATTGCTGCTCTGATTGCTTCCATGGGGATCTT  
 CTTGTTTTCAAGAACCTTAGCTGCCAGAGAGTACTCTGCACAAGAACATGTTCCCTCACTTATATCCTG  
 AACTCTATCATTATCATATCCACTGGTTGAGGTTGTGCCCAATGGAGACCTGGTGCGACGGGATCCTA  
 TGCATATTTTTATCATAACACATATATGTGGACAATGCAGTGGGA**ACTGTCACCACCCTTACCCCTGAG**  
 TGCACACGAGGAAAAGATGGACCCTCATGACAGTGAAGTGATAAGTTGCAAGATTCTGCACCTTTTCCAT  
 CAGTACATGATGGCTTGC**AACTACTTCTGGATGCTGTGTGAGGGGATCTATCTTCACACTCTGATCGTCA**  
 TGGCTGTGTTCA**CCGAGGATCAACGTCTGCGCTGGTACTATTTACTTGGCTGGGGTTCCCAATAGTGCC**  
 AACTATTATCCATGCCATTACTCGTGCCGTCTACTACAACGACA**ACTGTTGGCTGAGTACGGAGACCCAC**  
 TTGCTTTACATCATCCATGGACCCGTATGGCGGCTCTGGTGGTCAACTTCTTCTTCTGCTCAACATTG  
 TCCGTGTGCTTGTGACCAAGATGAGGCAAACCCATGAAGCCGAGGCCTACATGTACCTGAAGGCTGTAAA  
 GGCCACCATGGTCTTGTGCCCTGCTGGGGATCCAGTTTGTGTGTTTCCCTGGAGGCCCTCCAACAAA  
 GTGCTTGGGAAGATCTATGATTATCTCATGC**ACTCTGATTCA**TTTCCAGGGATTCTTTGT**CGCGACTA**  
 TCTACTGCTTCTGTAACCATGAGGTGCAAGTGACCCTGAAGCGTCA**GTGGGCGCAGTTCAAGATCCAGTG**  
 GAGCCATCGTGGGGTAGGCGCCGCCCCCAACCGAGTGGTTAGCGCTCCTCGTGTGTAGCCTTT  
 GCAGAACCGGTGGCTCCCATTTACATCTGCCATCAGGAACCAAGGAATCCTCCAGT**CAGCAACAACG**  
 AAGCGAAGAGGGTACTGAAATGATCCCATGAATGTCATCCAGCAAGACTCATCCGCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja1738\\_d06.zip](https://cdn.origene.com/chromatograms/ja1738_d06.zip)

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM\_053816

**Insert Size:** 1602 bp

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_053816.2</a> , <a href="#">NP_446268.2</a>
<b>RefSeq Size:</b>	3544 bp
<b>RefSeq ORF:</b>	1602 bp
<b>Locus ID:</b>	116506
<b>UniProt ID:</b>	<a href="#">P32214</a>
<b>Cytogenetics:</b>	4q13
<b>Gene Summary:</b>	<p>receptor for calcitonin; acts via adenylate cyclase signaling [RGD, Feb 2006] Transcript Variant: This variant (1b), previously called 1, represents the longer variant and encodes the longer isoform (1b, previously called a).</p>