

## Product datasheet for **RN200344**

### **Kcnip2 (NM\_020095) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kcnip2 (NM\_020095) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kcnip2  
**Synonyms:** Kchip2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200344 representing NM\_020095  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCGGGGCCAAGGCAGAAAGGAGAGTTTGTCCGAATCCCAGATCTGGACGGCTCCTATGACCAGCTTA  
CGGGCCACCTCCAGGGCCAGTAAAAAGCCCTGAAGCAGCGTTTCTCAAGCTGCTGCCGTGCTGCGG  
GCCCAAGCCCTGCCCTCAGTCAGTGAAAACAGCGTAGAGGATGAGTTTGAATTATCCACGGTGTGTAC  
CGACCTGAGGGCCTGGAACAACCTCCAGGAACAGACCAAGTTCACACGCAGAGAGCTGCAGGTCTGTACC  
GAGGCTTCAAGAACAATGCCCCAGTGGGATTGTCAACGAGGAGAACTCAAGCAGATTTATTCTCAGTT  
CTTTCCCAAGGAGACTCCAGCAACTATGCTACTTTTCTTTCAATGCCTTTGACACCAACCACGATGGC  
TCTGTCAGTTTTGAGGACTTTGTGGCTGGTTTGTCCGGTATTCTTCGGGGGACCATAGATGATAGACTGA  
GCTGGGCTTTCAACTTATATGACCTCAACAAGGACGGCTGTATCACAAGGAGGAAATGCTTGACATTAT  
GAAGTCCATCTATGACATGATGGCAAGTACACATACCCTGCCCTCCGGGAGGAGGCCCAAGAGAACAC  
GTGGAGAGCTTCTCCAGAAGATGGACAGGAACAAGGACGGCGTGGTGACCATCGAGGAATTCATCGAGT  
CTTGTC AACAGGACGAGAACATCATGAGGTCCATGCAGCTCTTTGATAATGTCATC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_020095  
**Insert Size:** 759 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_020095.2</a> , <a href="#">NP_064480.2</a>
<b>RefSeq Size:</b>	2316 bp
<b>RefSeq ORF:</b>	759 bp
<b>Locus ID:</b>	56817
<b>UniProt ID:</b>	<a href="#">Q9JIM59</a>
<b>Cytogenetics:</b>	1q54
<b>Gene Summary:</b>	Ca(2+)-binding protein that acts as an auxiliary subunit for voltage-gated K(+) channels [RGD, Feb 2006] Transcript Variant: This variant (b) lacks an exon in the 5' coding region compared to variant a. The resulting protein (isoform b) is shorter but has the same N- and C-termini compared to isoform a.