

## Product datasheet for **RN200362**

### Plaur (NM\_134352) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Plaur (NM\_134352) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Plaur  
**Synonyms:** Par; Plaur3; uPAR; uPAR-2; uPAR-3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN200362 representing NM\_134352  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCCTCTGAGGCGGCGGCTGCTGCTGGTGGTGGTGACGACCTGTGTCCAGCCTCCAGG  
GCCTGCGGTGCATACAGTGCAGAAAGTAACCAGGACTGCCTGGTGGAGGAGTGTGCCCTGGCCAGGACCT  
CTGCAGAACACGGTCTTCGGGAATGGGAAGACGCCGAGGAGCTGGAGGTGGTACCAGAGGCCTGTGC  
CACAAGGAGAAGACCAACAGGACCATGAGCTACCGCATGGGCTCGGTGATCGTCAGCCTCACCGAGACCG  
TGTGTGCCACAACTCTGCAACAGGCCAGACCCGGAGCCGAGGGCGACCTTTCCCGGGGCGGTTA  
CCTTGAGTGTGCATCCTGCACCTCTCTGGACCAGAGCTGTGAGAGGGGCCGGGAGCAAAGCTGCAGTGC  
CGGTATCCTACAGAACACTGTATTGAAGTGGTACGCTCCAAAGCACAGAACGGAGCGTGAAGGACGAAC  
CCTACACCAAAGGCTGTGGCAGTCTTCCCGGATGCCCGGGCACAGCAGGTTTCCATAGCAACCAGACCTT  
TCACTTCTGAAGTGTGCAACTTCACCAATGCAATGGTGGCCAGTCTGGATCTTACAGACTTACCA  
CCGAATGGCTTCCAATGTTACAGCTGTGAGGGGAACGAGCACCTTTGGATGTTCTACGAAGAGACGTCC  
TCATTGACTGCCGGGACCAATGAATCAGTCTTGGAGGCTACAGGCTTATAGTGTCTGGAAACCGGAG  
TTACACCGTAAGAGGCTGCGCCACGGCTTCTGGTGCCAAGGTTCCACGTGGCCGACTCCTTCCAGACC  
CAGTCAACCTCTCTATCTCTGCTGTAATGGCAGTGGCTGTAACCGCCCTACAGGGGGCGCCCCGGGC  
CAGGCCCTGCTCATCTTATCCTCATTGCCTCCCTGCTCCTGACCCCTCAGACTGTGGGCATCCCTCTCTG  
GACCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja1796\\_f06.zip](https://cdn.origene.com/chromatograms/ja1796_f06.zip)



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_134352
<b>Insert Size:</b>	987 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_134352.1</a> , <a href="#">NP_599179.1</a>
<b>RefSeq Size:</b>	1277 bp
<b>RefSeq ORF:</b>	987 bp
<b>Locus ID:</b>	50692
<b>UniProt ID:</b>	<a href="#">P49616</a>
<b>Cytogenetics:</b>	1q21
<b>Gene Summary:</b>	receptor for urokinase-plasminogen activator; may play a central role in malignant tumour progression [RGD, Feb 2006]