

## Product datasheet for **RN200490**

### **Prlr (NM\_012630) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Prlr (NM\_012630) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Prlr  
**Synonyms:** RATPRLR  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200490 representing NM\_012630  
**Red=Cloning site Blue=ORF Orange=Stop codon**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCCATCTGCACCTTGTCTTCGCTACTTGTCTCAACATCAGCCTCCTGAAGGGACAGTCACCACCAG  
GGAAACCTGAGATCCACAAATGTCGCTCCTGACAAGGAAACATTCACCTGCTGGTGAATCCTGGGAC  
AGATGGAGGACTTCTACCAATTATTCAGTACTTACAGCAAAGAAGGAGAGAAAACACCTACGAATGT  
CCAGACTACAAAACAGTGGCCCAACTCTGCTTCTTTAGCAAGCAGTACACTTCCATCTGGAAAATAT  
ATATCATCACAGTAAATGCCACGAACCAATGGGAAGCAGTTCCTCGGATCCACTTTATGTGGATGTGAC  
TTACATCGTTGAGCCAGAGCCTCCTCGGAACCTGACATTAGAAGTAAAACAGCTAAAAGACAAAAAACA  
TATCTGTGGTAAAAATGGTCCCCACCCACCATAACTGATGTGAAAACCTGGTTGGTTTACAATGGAATATG  
AAATTCGATTAAGCCTGAAGAAGCAGAAGAGTGGGAGATCCATTTTACAGGTCATCAAACACAGTTTAA  
AGTTTTTGACCTATATCCAGGGCAAAAGTATCTGTCCAGACTCGCTGCAAGCCAGACCATGGATACTGG  
AGTAGATGGAGCCAGGAGAGTTCGGTTGAAATGCCAAATGACTTCACCTTGAAGGACACAACCGTGTGGA  
TCATTGTGGCCATTCTCTGTGTCATCTGTTTATTGATTATGGTCTGGGCAGTGGCTTTGAAGGGCTATAG  
CATGATGACCTGCATCTTCCACAGTTCCTGGGCCAAAAATAAAAGGATTTGATACCCATCTGCTGGAG  
ACTGGTTCTCAAGCAAATACAAGGTGGATCTGTACCTGGCCCTTCTGGAGATTCCAGAACTGGACA  
ATGCAGGGGAACCTTGACTACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_012630  
**Insert Size:** 933 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_012630.2</a> , <a href="#">NP_036762.1</a>
<b>RefSeq Size:</b>	1633 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	24684
<b>UniProt ID:</b>	<a href="#">P05710</a>
<b>Cytogenetics:</b>	2q16
<b>Gene Summary:</b>	This is a receptor for the anterior pituitary hormone prolactin.[UniProtKB/Swiss-Prot Function]