

Product datasheet for **RN200495**

Pfkfb2 (NM_001033965) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pfkfb2 (NM_001033965) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pfkfb2
Synonyms:	PFK-2; RH2K
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >RN200495 representing NM_001033965
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGAGAATAGTACATTTTCCACAGAAGACAGCAGCAGCAGCTATAAACCCACGCTCAAACC
 TCCGAAGGGCAGGGAAGAAATGCTCATGGCTTCTACATGACCAACTCCCAACTCATTGTTATGAT
 TGGCTTCCAGCCCGGGTAAGACTTATGTGTCCAAGAAATTAACACGCTACCTGAACTGGATTGGAGTA
 CCTACCAAAGTGTAACTTCTGGTGTATATCGACGGGAAGCAGTCAAGTCTATAAGTCTATGACTTCT
 TTCGACATGACAATGAAGAGGCCATGAAGATCCGCAAACAGTGTGCCCTGGTGGCACTGGAAGATGTGAA
 GGCCTACTTTACTGAAGAGAGTGGCAGATCGCGGTGTTTGTGCCACCAATACCACTCGGGAGAGGAGG
 GACATGATTTTGAACCTTCCCAAGCAGAATGCCTTCAAGGTATTCTTTGTGGAATCTGTGTGTGATGATC
 CTGATGTCATTGCTGCCAATATTCTGGAGGTAAGGTGTCAAGCCCTGACTACCCCGAAAGGAATAGGGA
 GAATGTGATGGAGGACTTCTGAAGAGAATTGAGTGTACAAGGTCACTTACCAGCCCTTGACCCAGAC
 AACTATGATAAGGACCTCTCGTTCAAAAGGTGATGAATGTAGGCCAGAGGTTTCTGGTCAACAGAGTTC
 AGGACTACATCCAGAGTAAGATTGTCTACTACCTGATGAACATCCATGTCCATCCTCGCACCATCTATCT
 GTGCCGGCACGGAGAGCGAATTAATCTTTGGGAAAGATTGGGGGTGACTCTGGCCTTTCTGTTGCGA
 GGAAAGCAGTTTCTCAGGCTCTGAAGAAGTTTCTGGAGGAACAGGAGATCCAGGACCTCAAAGTGTGGA
 CAAGCCAGTTGAAGAGGACAATTCAGACTGCTGAGTCTCTTGGGGTGACCTATGAGCAATGGAAGATCCT
 AAATGAGATTGATGCTGGCGTGTGTGAGGAGATGACTTATTCGGAGATCGAACACGGTATCCAGAGGAA
 TTTGCACTTCGAGATCAAGAGAAGTATCTGTATCGATATCCTGGTGGGGAGTCATACCAGGACCTGGTGC
 AGCGGCTGGAGCCTGTGATCATGGAGCTGGAGCGGCAAGGCAACGTCCTCGTTATCTCACCAGGCTGT
 CATGCGCTGCCTCTGGCCTACTTCTGGATAAAGGCGCAGATGAGTTGCCGTACTTGAGGTGCCCCCTC
 CACATCATCTTCAAACCTACTCCTGTGGCCTATGGTTGCAAAGTGGAAACAATTACACTGAATGTGGAAG
 CTGTGGACACACATCGTGACAAGCCAACCTCACAACCTCCCAAGAGCCAAACCCCTGTAAGGATGAGAAG
 GAACAGCTTTACGCCTCTGTCCAGTTCGAACACAATAAGGCGTCCAAGAAATTACAGTGTGGGAGCCGA
 CCCCTCAAGCCCTCAGCCCTCTCCGTGCCCTGGATATGCAAGAAGGGGCGGACCAGCCGAAGACCCAAG
 TCCAGCAAGGGTCAGCCAGGCTACAGAACATCTACAGAAAGCCCTGGAGTTTCAAATGGACATAGGGA
 AGTGGAGAATGTGCTGGTAAGCATAGACGCCCTCAATGGCATCCCTCACTTTGCTCTCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** Sgfl-Mlul
- ACCN:** NM_001033965
- Insert Size:** 1674 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033965.2](#), [NP_001029137.1](#)

RefSeq Size: 3824 bp

RefSeq ORF: 1674 bp

Locus ID: 24640

UniProt ID: [Q9JH5](#)

Cytogenetics: 13q13

Gene Summary: catalyzes the synthesis and degradation of fructose 2,6-bisphosphate [RGD, Feb 2006]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.