

Product datasheet for **RN200530**

Acp1 (NM_021262) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Acp1 (NM_021262) Rat Untagged Clone
Tag: Tag Free
Symbol: Acp1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN200530 representing NM_021262
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGAGGTTGGGTCCAAGTCAGTGCTGTTCTCGTGTCTCGGTAACATTTGCCGGTCACCCATTGCAG
AAGCAGTGTTTCAGAAAATTGGTAACTGATGAAAACGTTTCAGATAACTGGAGGATAGACAGTCCGGCTAC
ATCCACCTATGAAGTGGGAACCCCTGACTATCGAGGGCAGAAGTGCATGAAAAACATGGCATCCAC
ATGCAACACATTGCACGGCAGATTACAAGAGAAGACTTTGCCACATTTCGATTATATACTGTGTATGGATG
AAAGCAATCTGAGAGATCTGAATAGAAAAAGTAATCAAGTAAAAACTGCAAAGCTAAAAATCGAGCTACT
TGGGAGCTATGATCCACAGAAGCAGCTCATTATTGAAGATCCCTATTATGGCAATGACTCTGACTTCGAG
GTGGTGTACCAGCAATGCCTTAGGTGCTGCAAGGCCTTCTGGAGAAGACTCACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_021262
Insert Size: 477 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021262.2](#), [NP_067085.1](#)

RefSeq Size: 1436 bp

RefSeq ORF: 477 bp

Locus ID: 24161

UniProt ID: [P41498](#)

Cytogenetics: 6q16

Gene Summary: catalyzes the hydrolysis of p-nitrophenyl phosphate; may play a role in synaptic transmission [RGD, Feb 2006]

Transcript Variant: This variant (1) represents the shorter transcript.