

Product datasheet for **RN200567**

Hmg20b (NM_001108731) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hmg20b (NM_001108731) Rat Untagged Clone
Tag: Tag Free
Symbol: Hmg20b
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN200567 representing NM_001108731
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCACGGTCCCCGGCAGCCAGGGCAGCCACCGCACCTGCGGGAGGCAAGACTCCGGGCCAGCATG
 GGGCCTTTGTAGTGGCTGTCAAGCAGGAGCGCAGCGAAGGTCCCCGTGCCGAGAGAAGGGGCCAGGA
 GGAGGAGCCAGTGAAGAAGCGAGGCTGGCCAAAGGCAAGAAAAGAAAGAAGATTCTCCCAATGGGCC
 AAGGCCCGGTACGGGCTATGTTTCGCTTCCTGAACGAGAGGGCGGAGCAGATTCGCACCCGACATCCGG
 ACCTTCCCTTCCCGGAGATCACCAAGATGTTGGCGCGGAGTGGAGCAAGCTGCAGCCAGCAGAGAAGCA
 GCGGTACCTGGATGAGGCAGAGAAGGAGAAACAGCAGTACTTGAAGGAGCTGTGGGCATACCAGCAGTCA
 GAGGCCTACAAGTCTGCACCGAGAAGATCCAAGAAAACAAAATAAAAAAGAGGACTCCAGCTCTGGGC
 TCATGAACACCCTCCTGAATGGACACAAGGGTGTGGACTGCGGCGATGGCTTTTCAACCTTTGATGTTCC
 CATCTTCACTGAAGAGTTTTTGGACAAAACAAAGCACGCGAGGCAGAGCTGCGGCGCTGCGGAAGATG
 AACGTGGCTTTCGAGGAGCAGAACGCAGTGTGCAGCGCCACACTCAGAGCATGAACAGCGCTCGTGAGC
 GTCTGGAGCAAGAGCTGGCGCTGGAGGAGCGGCGCAGCTGGCGTGCAACAGCAACTGCAGGCTGTGAG
 ACAAGCGCTCACGGCTAGTTTCGCTCTCTGCCGGTCCCCGGCACCGGGAGACGCCGACGCTCGGAACG
 CTGGACTTCTACATGGCGGGCTTACGGAGCCATCGAGCGGACCCGGCGCAGCAGAGAGACTGATCG
 CGCGCGTCAAGGAGATCCTGGCGGAGTAGCCAGCGAGCACCT**GTGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001108731
Insert Size: 957 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108731.1</u> , <u>NP_001102201.1</u>
RefSeq Size:	1676 bp
RefSeq ORF:	957 bp
Locus ID:	362825
Cytogenetics:	7q11