

Product datasheet for **RN200593**

Sh3glb2 (NM_001009692) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sh3glb2 (NM_001009692) Rat Untagged Clone
Tag: Tag Free
Symbol: Sh3glb2
Synonyms: MGC105723
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN200593 representing NM_001009692
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGGAGCTGCTGGCTCAGTACATGGCAGAGGCAGCCAGCGAACTGGGCCCAACTCCCTATGAGTG
 GCTTGTGCTCTGCCCTCTCACCTATGACTCTTCCAGGGAAGACTGATGAAGGTGTCAGAAGC
 TGAAGCGTCTGGGAGCCGAGAACGGGATTTTCATTCACTGCCTCCCTCAGTTCTCACACCCCTTG
 CGCAACTTCTGGAAGGGGACTGGAAAACGATTCGAAGGAGAGGCGGCTTCTGCAGAACCAGCGCCTTG
 ACTTGGATGCCTGCAAAGCCCGCTAAAGAAGGCCAAGGCGCAGAAAGCCAAAGCCACGTGTGAGGGAGA
 TACGGTGCCTGACTTTCAGGAGACTAGACCTCGTAATTACATTCTCTCGGCCAGCGCCTCCGCGCTTTGG
 AACGATGAGGTAGACAAGGCTGAGCAGGAGCTTCGAGCGGCCAGACCGAGTTTGACCGGCAGGCAGAAG
 TTAAGTCTCAGACCACCTACTATGCACAGTGCTACCGCCACATGCTGGATCTGCAGAAACAGCTAGGCAGC
 TCCAGGGTGCCATATTTCCAGGCACCTTCGTGGGCACTACAGAGCCTGCCTCCCACCCCTGAGCAGCA
 CCTCACCCACCACCACTGCGGCCACCATGCCTGTGGTACCTACTGGGGCTGGCTTGGCCCTCCAGAAGA
 GGCAGCACTGTGCTGGAGGAGGTGGCTCCACCAGCGGGACTCGAAAGGCCCGGGTGTCTACGAC
 TATGAGGCAGCTGACAGCAGTGAGCTGGCCCTGTGGCTGACGAGCTCATCACTGTCTACAGCCTGCCGG
 GCATGGACCCCGACTGGCTCATTGGTGAGAGGCAACAAGAAGGGCAAGGTTCTGTCCCTACCTGGA
 ACTTCTCAGTTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001009692



[View online »](#)

Insert Size:	993 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001009692.1</u> , <u>NP_001009692.1</u>
RefSeq Size:	1849 bp
RefSeq ORF:	993 bp
Locus ID:	311848
UniProt ID:	<u>Q5PPI9</u>
Cytogenetics:	3p12