

## Product datasheet for **RN200631**

### Twist1 (NM\_053530) Rat Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Twist1 (NM\_053530) Rat Untagged Clone  
 Tag: Tag Free  
 Symbol: Twist1  
 Synonyms: Twist  
 Mammalian Cell Selection: Neomycin  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 Fully Sequenced ORF: >RN200631 representing NM\_053530  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGATGCAGGACGTGTCCAGCTCGCCAGTCTCGCCGGCCGACGACAGCCTGAGCAACAGCGAGGAGGAGC  
 CGGACCGGCAGCAGCCGGCGAGCGGCAAGCGGGGGCTCGCAAGAGACGCAGCAGTCGGCGCAGCGCGGG  
 CGGCAGCGCGGGCCCGCGGGCCACGGGCGGGGCATCGGAGGCGGCAGGAGCCGGGCAGCCCGGCC  
 CAGGGCAAGCGCGCAAGAAATCTCGGGCGGAGGCGGGCGCGGGCGGAGGCGGGCGGGCGGGCGGCA  
 GCAGCAGCGGGGGCGGGAGCCCGCAGTCGTACGAGGAGCTGCAGACACAGCGGGTCATGGCCAACGTGCG  
 GGAGCGCCAGCGACCCAGTCGCTGAACGAGGCATTTGCCGCCCTGCGCAAGATCATCCCCACGCTGCC  
 TCGGACAAGCTGAGCAAGATTGAGACCCTCAAAGTGGCGGCCAGGTACATCGACTTCTGTACCAGGTCC  
 TGCAGAGCGACGAGCTGGACTCCAAGATGGCAAGCTGCAGCTATGTGGCCACGAGCGGCTCAGTACGC  
 CTTCTCCGTCTGGAGGATGGAGGGGGCCTGGTCCATGTCCGCGTCCCAC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul  
 ACCN: NM\_053530  
 Insert Size: 612 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_053530.2](#), [NP\\_445982.1](#)

**RefSeq Size:** 1513 bp

**RefSeq ORF:** 612 bp

**Locus ID:** 85489

**Cytogenetics:** 6q16

**Gene Summary:** beta helix loop helix transcription factor [RGD, Feb 2006]