

Product datasheet for **RN200638**

Fgf17 (NM_019198) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fgf17 (NM_019198) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fgf17
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN200638 representing NM_019198 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGGAGCCGCCGCTGCTGCCTAACCTTACCCTGTGCTTGCAGCTATTGATCCTCTGCTGTCAAACGC
AGGGGGAGAATCACCCGCTCTCCTAATTTAAACCAGTACGTGAGGGACCAGGGCGCTATGACTGACCAGCT
GAGCAGGCGGCAAATCCGTGAGTACCAGCTCTACAGCCGACCAGTGCCAAGCACGTGCAGTCCACGGG
CGTCGCATCTCTGCCACCGCTGAGGATGGCAACAAGTTCGCCAAGCTCATCGTGGAGACAGACACATTCCG
GCAGCAGAGTCCGCATCAAGGGGCAGAGAGCGAGAAGTATATCTGTATGAACAAGAGGGCAAGCTGAT
TGGGAAGCCGAGCGGGAAGAGCAAAGACTGCGTATTCACCGAGATCGTACTGGAGAACAACACACAGCC
TTCCAGAACGCCCGCATGAGGGCTGGTTCATGGCTTTCACCCGGCAGGGCCGGCCTCGCCAGGCCCTCC
GGAGCCGCCAGAACCAGCGAGAGGCCACTTCATCAAGCGCTCTACCAAGGCCAGCTGCCTTCCCAA
CCACGCTGAAAGGCAGAAGCAGTTCGAATTCGTGGGCTCCGCCCCACTCGCAGGACCAAGCGACTCGG
AGGCCCCAGTCCCAAACGT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_019198
Insert Size:	651 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019198.1</u> , <u>NP_062071.1</u>
RefSeq Size:	651 bp
RefSeq ORF:	651 bp
Locus ID:	29368
UniProt ID:	<u>P63076</u>
Cytogenetics:	15p11
Gene Summary:	may mediate induction and patterning of the embryonic brain [RGD, Feb 2006]