

Product datasheet for **RN200651**

Cfap77 (NM_001127563) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Cfap77 (NM_001127563) Rat Untagged Clone

Tag: Tag Free

Symbol: Cfap77

Synonyms: RGD1564114

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN200651 representing NM_001127563
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGCTGACCCGGGGAAACCCCGAAAGACCTCTCGCATGGAAGAAGAAGAAGCAGCCCGTGACCCGCA
 CGGTGAGCCAGATCTGCCCGCCGCCCGACGGCCGCTGACGGTGGTCGACATCCGCCGGGCATGGAGAA
 CGAAAGGCTGGGGTCTGCGTGACTCCATGTTTCAGAACCCGCTCATTGTCAAGGCAGAACTTGGCAAG
 CCCAGGGAGAGAAGCTGCAGCCTGCCTGGCGTCGACTTTAATTATGGACTCTACATCCAAGGGATAGATG
 GAGGGTTCTGAAGCCATTGGGCACTGGAATGTGTTTAAGCAGCAACCCACCTGTCCCATGAGCTGAC
 CCGGAATATATTGCCATGAACCGAGGGGCACTCAAAGCTGGCCTGGTGACTGCCCGGGAAACATGCTC
 TATCGTGAATGAATGACATTGCAATCAATGACCAGGAAGACCGGCGACAGAAGGAGCCGCCCTCCATCC
 CTCCCAACATGACCTTTGGCATCCGTTACGGCCCTCAACACCCTTCTTTGACCTGTTACAACATCGGTA
 CCAGCAGCTATGGATACAGGAGCAGAAGGCCACTCAGCAAGCCATCAAAATGGAGAAGAAACAGAAGGTA
 ATCCTGGGAAAGCTGTATGAGACGAGGAGCAGTCAGCTGAGGAAGTACAAGCCTCCTGTGAAGCTGGATA
 CTCTCTGGCACATGCCTCACTTCAAGAAGGCGGTTTCCCATCTTGATACGTTCCCTACGGATGCCGCTCG
 CCAGAGGGCTCTCAAAGCCACAAGGAAGAGCATGCTGTACGCCAGGGCACCTACGAATGGGCAACTAC
 ACCCACCCT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001127563

Insert Size: 852 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001127563.1, NP_001121035.1</u>
RefSeq Size:	1739 bp
RefSeq ORF:	852 bp
Locus ID:	499765
Cytogenetics:	3p12