

Product datasheet for **RN200706**

Tp53inp1 (NM_181084) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tp53inp1 (NM_181084) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tp53inp1
Synonyms:	Stinp; Trp53inp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN200706 representing NM_181084 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTTCAGAGACTGAATAAAATGTTTGTAGGTGAAGTCACTACTTCTTCCAGCCACGAGCCAGAATTTA
 GTGAGAAAGAAGATGATGAGTGGATTCTTGTTGACTTCATAGATACCTGCACTGGTTTCTCAGCAGAAGA
 AGAGGAGGAAGATGAGGACATTGGTGAAGAGTCCTCAGCAGAGCACACTTCTGTCTTCTCCTGTTTACCT
 ACATCTTTGGAGTGCTTGGCTGACACCAGCGATTCTGCTTCTTCAGTTTGAGTCCTGTCCAATGGAGG
 AGAGCTGGTTCATCACTCTCCCCATGTTTACAGCAGGGGGATTAACCACTATCAAAGTGAAACAAG
 TCCTATGGAAAACCTTCTCATTGAACATCCCAGCATGTCTGTCTATGCTGTGCACAACTCTGTCTGGT
 CTCAGTGAAGCTAGCTGTGGGAACGATGACTATAACTCAAGTGGTCCCAGAATGGAAGCCAGAGTGAAC
 TGGGGACGCACATTCACTGCTGTGTGTCAGCGCTTGCTGCTCAGGCAACCTTTTGGAAACAACCCAGAG
 TTTTCGCCCTTCCCAGTGGATAAAAGGACACAGTGAAAGACAGTCTCTGAACAGAAATGGCTCCGTCGC
 CAGAACCTTACCAGGGATTGCCACACTCGACAAATGAAGCACAGGTGCTGGGCGGTTACCAGCCCTGCC
 CACGCCAGTACAATTAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_181084
Insert Size:	720 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_181084.2</u> , <u>NP_851598.1</u>
RefSeq Size:	1013 bp
RefSeq ORF:	720 bp
Locus ID:	297822
UniProt ID:	<u>Q80YE2</u>
Cytogenetics:	5q13
Gene Summary:	displays increased expression in response to stress [RGD, Feb 2006]