

## Product datasheet for **RN200770**

### **Cd82 (NM\_031797) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cd82 (NM\_031797) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Cd82  
**Synonyms:** Kai1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200770 representing NM\_031797  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGGGCAGGCTGTGTCAAAGTCACCAAGTACTTCTCTTCTTCAATCTGCTGTTCTTTATCCTGG  
GTGCTGTGATCCTGGGCTTTGGGGTGTGGATTCTGCAGACAAGAGCAGCTTCATTTCTGCCTACAAC  
CTCATCCAGCTCGCTGCAGGTGGGGGCATACGTCTTCATTGGTGTGGGCGCCATCACCATGCTGATGGGC  
TTCTGGGCTGTATTGGTGTGTCAATGAGGTCCGCTGCCTGCTGGGTCTGTACTTCGTCTTCTTCTCTGC  
TGATTCTTATAGCACAGGTGACTGTAGGGGTCCTTCTACTTCAATGCTGACAAGCTGAAGCAGGAGAT  
GGGGAACACAGTGATGGACATTATTCAGAACTACAGCGTCAATGCCAGCAGCAGCCGCGAGGAGGCTTGG  
GACTATGTGCAGGCGCAGGTGAAGTGTGTGGCTGGGTGAGCCCTCCAAGTGGACACGTAACCCGGTGC  
TCAAGAACTCTACCAAGACCACTTACCCATGCTCCTGTGAGAAGACGAAGGAGGAGGACAACCAGCTCAT  
TGTGAAGAAAGGGTTCTGTGAGTCTGACAACAGCACCCGAAGTAAAAACAGTCTGAGGACTGGCCTGTG  
CACCCGGAGGGGTGCATGGAGAAGGCTCAGGCGTGGCTACAGGAGAACTCGGCATCCTTCTGGGCGTGT  
GTGCTGGTGTGCTGCTATTGAGCTGCTGGGGTGTTCCTGTCCATATGTTTGTGCCGGTACATTCATTC  
TGAAGATTACAGCAAGGTCCCAAGTACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_031797  
**Insert Size:** 801 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_031797.2</a> , <a href="#">NP_113985.2</a>
<b>RefSeq Size:</b>	1838 bp
<b>RefSeq ORF:</b>	801 bp
<b>Locus ID:</b>	83628
<b>Cytogenetics:</b>	3q31
<b>Gene Summary:</b>	a metastasis suppressor gene; loss of expression is associated with increased metastatic potential of prostate cancer cells [RGD, Feb 2006]