

## Product datasheet for RN200776

### Hnrnph3 (NM\_001108532) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hnrnph3 (NM\_001108532) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Hnrnph3  
**Synonyms:** Hnrph3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN200776 representing NM\_001108532  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGC**C

ATGGATTGGGTTATGAAACATAATGGTCCAATGACGCTAGTGATGGGACAGTGCAGCTTCGTGGACTGC  
 CATTGGTTGCAGCAAAGAGGAAATAGTTCAGTCTTTCAAGGGTTGAAATCGTGCCAAATGGGATAAC  
 ATTGACGATGGACTACCAGGGGAGAAGCACAGGGGAGGCTTCGTGCAGTTTGCTCAAAGGAGATAGCA  
 GAAATGCTCTGGGAAACACAAGGAAAGAATAGGGCACAGGTATATTGAGATCTTCAGGAGTAGCAGGA  
 GTGAAATCAAAGGATTTATGATCCACCAAGAAGATTGCTGGGACAGCGACCTGGACCATATGATAGACC  
 AATAGGAGGAAGAGGGGGTTATTATGGAGCTGGGCGTGGAAAGTATGTATGACAGAATGCGACGAGGAGGT  
 GATGGATATGATGGTGGCTATGGAGTTTTGATGACTATGGTGGCTATAATAATTACGGCTATGGGAATG  
 ATGGCTTTGATGACAGGATGAGAGATGGAAGAGGCATGGGAGGACATGGCTATGGTGGAGCTGGTATGC  
 AAGTTCAGGTTCCATGGTGGTCATTTGTGCACATGAGAGGGCTGCCTTTTCGTGCGACTGAAAATGAT  
 ATTGCTAATTTCTTCACCACTCAATCCAATACGCGTTCATATTGATATTGGAGCCGATGGCAGAGCAA  
 CAGGAGAGGCAGATGTAGAGTTGTGACACATGAAGATGCAGTGGCTGCCATGTCTAAAGATAAAAACAA  
 CATGCAACATCGATACATTGAACTCTTCTTGAATTCACCCCTGGAGGTGGCTCTGGAATGGGAGGTTCT  
 GGAATGGGAGGCTATGGCAGAGATGGAATGGATAATCAAGGAGGCTACGGATCTGTTGGGAGGATGGGAA  
 TGGGGAACAGCTACAGTGGAGGCTATGGTCTCCTGATGGCTGGGAGGCTATGGTCTGGTGGTGGAGG  
 CAGTGGAGGTTACTATGGCAAGGTGGAATGAGTGGAGGTGGATGGCGTGGGATGTACT**GA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-RsrII  
**ACCN:** NM\_001108532



[View online »](#)

<b>Insert Size:</b>	1041 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001108532.1</a></u> , <u><a href="#">NP_001102002.1</a></u>
<b>RefSeq Size:</b>	1879 bp
<b>RefSeq ORF:</b>	1041 bp
<b>Locus ID:</b>	361838
<b>Cytogenetics:</b>	20q11